

Thermo Scientific Conventional Restriction Enzymes Buffer Activity Chart

Enzyme	Specificity 5'→3'	Recommended buffer for 100% activity	Enzyme activity in Thermo Scientific™ buffers, %					Thermo Scientific™ Tango™ buffer for double digestion	Enzyme properties	
			B (blue) 1X	G (green) 1X	O (orange) 1X	R (red) 1X	Tango™ (yellow) 1X			
AanI (PstI)	TTA_TAA	Tango	50-100	50-100	0-20	0-20	100	20-50	1X or 2X	☺
AarI	CACCTGC(4/8)↓	Unique	NR	NR	0-20	0-20	NR	50-100	2X	☺
AatII	GACGT↓C	Tango	50-100	20-50	0-20	0-20	100	20-50	1X or 2X	☺
Acc65I (Asp718I)	G1GTACC	O	0-20	20-50	100	20-50	20-50	50-100	1X or 2X	☺
AdelI (DraII)	CACNNN↓GTG	G	0-20	100	20-50	100	100*	20-50	1X* or 2X	☺
AjilI (BmgBI)	CAC↓GTC	Unique	NR	NR	20-50*	NR	NR	20-50*	2X*	☺
AjuI	↓(7/12)G(A)N,TTGG(11/6)↓	R	0-20	50-100	20-50	100	50-100	50-100	1X or 2X	☺
AluI	↓(10/12)G(A)N,TCG(12/10)↓	R	0-20	50-100	20-50	100	50-100	50-100	1X or 2X	☺
AluI (30°C)	↓(7/12)G(A)N,TCG(12/10)↓	R	0-20	50-100	20-50	100	20-50	100	1X or 2X	☺
AluI	AG↓CT	Tango	50-100	0-20	0-20	0-20	100	20-50	1X or 2X	☺
Alw211 (BstHKA1)	SWGCC↓C	O	0-20	20-50	100	50-100	20-50	50-100	1X or 2X	☺
Alw261 (BsmAI)	GTCTC(1/5)↓	Tango	50-100	100	0-20	0-20	100	100	1X or 2X	☺
Alw441 (ApaLI)	G1GTACC	Tango	50-100	100	0-20	50-100	100	50-100	1X or 2X	☺
ApaI	GGGCC↓C	B	100	20-50	0-20	0-20	20-50	0-20	1X	☺
BamHI	G↓GATCC	Unique	20-50*	100	20-50	50-100*	100*	50-100	1X* or 2X	☺
BauI (BssSI)	CACGAG↓(5/1)↓	Tango	0-20	50-100	0-20	50-100	100	50-100	1X or 2X	☺
BclI (55°C)	T↓GATCA	G	20-50	100	20-50	20-50	100*	100	1X* or 2X	☺
BclI (NotI)	C↓GGG	Tango	20-50	50-100	50-100	50-100	100	50-100	1X or 2X	☺
BclI (SpeI)	A↓CTAGT	Tango	50-100	50-100	0-20	20-50	100	0-20	1X	☺
BfmI (SfiC)	C↓TFRAG	Tango	0-20	50-100	0-20	0-20	100	20-50	1X or 2X	☺
BfuI (BciVI)	GTATCC(6/5)↓	Unique	NR	NR	0-20	0-20	NR	50-100*	2X*	☺
BglI	GCNNNN↓NGGC	O	0-20	50-100	100	100	0-20	100	2X	☺
BglII	A↓GATCT	O	0-20	20-50	100	50-100	0-20	100	2X	☺
BmeI390I (SrfAI)	CC↓NGG	O	20-50	50-100	100	50-100	50-100	50-100	1X or 2X	☺
BoxI (PstAI)	GACNN↓NGTCC	Tango	0-20	0-20	0-20	20-50	100	20-50	1X or 2X	☺
BpiI (BbsI)	GAGAC(2/6)↓	G	20-50	100	50-100	50-100	50-100	50-100	1X or 2X	☺
BplI	↓(8/13)G(A)N,CTC(13/8)↓	Tango	0-20	20-50	0-20	0-20	100	20-50	1X	☺
Bpu10I	CCTNAGC(1-5/2)↓	Unique	0-20	20-50*	50-100*	100*	50-100*	100*	1X* or 2X*	☺
Bpu102I (BplI)	GC↓TAGC	Tango	50-100	50-100	20-50	20-50	100	20-50	1X or 2X	☺
BseDI (BsaJI) (55°C)	C↓CNNNGG	Tango	50-100	20-50	0-20	0-20	100	50-100	1X or 2X	☺
BseGI (BtsCI) (55°C)	GGATG(2/0)↓	O	20-50	50-100	20-50	20-50	100	20-50	1X or 2X	☺
BseJI (BsaBI) (65°C)	GATNN↓NATCC	O	NR	100*	0-20	NR	NR	100*	2X*	☺
BseLI (BsaCI) (55°C)	CCNNNN↓NNGG	Tango	20-50	100	50-100	20-50	100	50-100	1X or 2X	☺
BseMI (BspCNI) (55°C)	GCATG(2/0)↓	R	0-20	20-50	0-20	100	50-100	50-100	1X or 2X	☺
BseNI (BsrI) (65°C)	CTCAG(10/8)↓	Tango	50-100	50-100	50-100	50-100	50-100	50-100	1X or 2X	☺
BsePI (BspCNI) (55°C)	CTCAG(10/8)↓	Tango	50-100	50-100	50-100	50-100	50-100	50-100	1X or 2X	☺
BseXI (BmeI580I) (55°C)	GKGCAC↓C	G	20-50	100	0-20	20-50	50-100	0-20	1X	☺
BseXI (BbvI) (65°C)	GCAGC(8/12)↓	Unique	NR	NR	NR	NR	NR	NR	NR	☺
Bsh1236I (BstUI)	CG↓CG	R	0-20	0-20	50-100	100	20-50	50-100	1X or 2X	☺
Bsh1285I (BseI)	CGR↓CG	G	20-50	100	20-50	0-20	0-20	20-50	2X	☺
BshNI (BamI)	G↓GATCC	O	0-20	20-50	100	50-100	0-20	100	2X	☺
BshTI (AgeI)	A↓CCGGT	O	0-20	20-50	100	50-100	20-50	20-50	1X or 2X	☺
Bsp68I (NruI)	TCG↓CGA	O	0-20	20-50	100	50-100	20-50	50-100	1X or 2X	☺
Bsp119I (BstBI)	TT↓CGAA	Tango	20-50	0-20	0-20	0-20	100	100	1X or 2X	☺
Bsp120I (PvuOMI)	G↓GGCC	B	100	20-50	0-20	20-50	50-100	0-20	1X	☺
Bsp143I (Sau3AI)	↓GATC	Unique	20-50	20-50	0-20	0-20	50-100	20-50	1X or 2X	☺
Bsp1407I (BsrGI)	T↓GTACA	Tango	0-20	20-50	0-20	20-50	100	50-100	1X or 2X	☺
BspLI (NlaIV)	GGN↓NCC	Tango	50-100	50-100	0-20	20-50	100	20-50	1X or 2X	☺
BspOI (BmtI)	GCTAG↓C	O	0-20	0-20	100	100	0-20	20-50	2X	☺
BspP1 (AlwI) (55°C)	GGATC(4/2)↓	Tango	20-50	20-50	0-20	0-20	100	0-20	1X	☺
BspT1 (AflII)	C↓TTAAG	O	0-20	0-20	100	20-50	0-20	50-100	2X	☺
Bst1107I (Bst217I)	GTA↓TAC	O	20-50	50-100	100	100	20-50	100	1X or 2X	☺
BstXI (55°C)	CCNNNN↓NTGG	O	0-20	50-100	100	50-100	100	100	1X or 2X	☺
Bsu15I (ClaI)	AT↓CGAT	Tango	20-50	20-50	20-50	20-50	50-100	20-50	1X or 2X	☺
BsuRI (HaeIII)	GG↓CC	R	20-50	20-50	50-100	100	50-100	100	1X or 2X	☺
BveI (BspMI)	ACCTGC(4/8)↓	O	0-20	20-50	100	20-50	50-100	100	1X or 2X	☺
CalI (AlwNI)	CAGNN↓CTG	Tango	20-50	20-50	20-50	50-100	100	50-100	1X or 2X	☺
Cfr9I (XmaI)	C↓CCGGG	Unique	0-20	0-20	0-20	0-20	20-50	0-20	1X	☺
Cfr10I (BsrFI)	R↓CCGGY	Unique	0-20	20-50	20-50	50-100*	20-50	50-100	1X or 2X	☺
Cfr13I (Sau36I)	G↓GNCC	Tango	50-100	50-100	20-50	20-50	100	20-50	1X or 2X	☺
Cfr42I (SacI)	CCGC↓GG	B	100	50-100	0-20	0-20	50-100	0-20	1X	☺
CpoI (RsrII)	CG↓GWCCG	Tango	20-50	50-100	50-100	20-50	100	50-100	1X or 2X	☺
CseI (HgaI)	GAGCC(5/10)↓	R	NR	50-100*	50-100	100*	50-100	1X* or 2X*	☺	
Csp6I (CviOI)	G↓TAC	B	100	50-100	0-20	0-20	50-100	0-20	1X	☺
DpnI	GmE↓ATC	Tango	100	100	50-100	50-100	100	50-100	1X or 2X	☺
DraI	TTT↓AAA	Tango	50-100	50-100	20-50	20-50	100	50-100	1X or 2X	☺
Eam1104I (EarI)	CTCTT(1/4)↓	Tango	50-100	50-100	0-20	0-20	100	0-20	1X	☺
Eam1105I (AhdI)	GACNN↓NNGTCC	Unique	20-50	50-100	0-20	0-20	50-100	20-50	1X or 2X	☺
Ecl136I (EcoRI)	GAG↓CTC	Unique	50-100	20-50	0-20	0-20	50-100	0-20	1X	☺
Eco24I (BamI)	GRG↓C	Tango	50-100	50-100	0-20	20-50	100	0-20	1X	☺
Eco31I (BsaI)	GGTCTC(1/5)↓	G	50-100	100	0-20	0-20	50-100	20-50	1X or 2X	☺
Eco32I (EcoRV)	GAT↓ATC	R	0-20	50-100	50-100	100	20-50	100	1X or 2X	☺
Eco47I (AflI)	G↓GWCC	R	0-20	50-100	50-100	100	50-100	50-100	1X or 2X	☺
Eco47II (AvaI)	AG↓GCT	O	0-20	20-50	100	100	50-100	100	1X or 2X	☺
Eco52I (EagI)	C↓GGCCG	Unique	0-20	0-20	0-20	20-50	0-20	20-50	2X	☺
Eco57I (AclI)	CTGAAG(16/14)↓	G	100	100	20-50	20-50	50-100	1X or 2X	☺	
Eco72I (PmlI)	CAC↓GTC	Tango	NR	NR	0-20	0-20	100	20-50	1X or 2X	☺
Eco81I (Bsu36I)	CC↓TAGG	Tango	50-100	100	0-20	0-20	100	0-20	1X	☺
Eco88I (AvaI)	C↓YCGRG	Tango	100	50-100	0-20	0-20	100	20-50	1X or 2X	☺
Eco91I (BstEII)	G↓GTNACC	O	20-50	20-50	100	50-100	NR	100	1X or 2X	☺
Eco105I (SnaBI)	TAC↓GTA	Tango	100*	50-100	0-20	0-20	100	0-20	1X	☺
Eco130I (SfiI)	C↓CWGGG	O	0-20	20-50	100	50-100	50-100	100	1X or 2X	☺
Eco147I (SfiI)	AGG↓CCT	B	100	50-100	20-50	20-50	50-100	0-20	1X	☺
Eco109I (DraI)	RG↓NCCY	Tango	50-100	20-50	20-50	20-50	100	100	1X or 2X	☺
EcoRI	G↓AATTC	Unique	0-20	NR	100	100*	NR	100	2X	☺
EcoRII	↓CCWGG	O	20-50	50-100	100	50-100	20-50	50-100	1X or 2X	☺
EheI (SbfI)	GGC↓GCC	Tango	20-50	50-100	0-20	0-20	100	20-50	1X or 2X	☺
Esp3I (BsmBI)	CGTCTC(1/5)↓	Tango	100	20-50	0-20	0-20	100	0-20	1X	☺
FaqI (BsmFI)	GGGAC(10/14)↓	Tango	20-50	20-50	0-20	0-20	100	20-50	1X or 2X	☺
FspAI	RTGC↓GAY	O	0-20	0-20	100	50-100	0-20	50-100	2X	☺
FspBI (SfaI)	C↓TAG	Tango	50-100	20-50	0-20	0-20	100	0-20	1X	☺
GsuI (BpmI) (30°C)	CTGGAG(16/14)↓	B	100	50-100	20-50	20-50	100	50-100	1X or 2X	☺
HhaI	GCG↓C	Tango	50-100	50-100	20-50	20-50	100	20-50	1X or 2X	☺
Hin1I (BsaHI)	GR↓CGYC	G	20-50	100	20-50	20-50	20-50	20-50	1X or 2X	☺

Enzyme	Specificity 5'→3'	Recommended buffer for 100% activity	Enzyme activity in Thermo Scientific™ buffers, %					Thermo Scientific™ Tango™ buffer for double digestion	Enzyme properties	
			B (blue) 1X	G (green) 1X	O (orange) 1X	R (red) 1X	Tango™ (yellow) 1X			
Hin1II (NlaII)	CATG↓	G	50-100	100	20-50	50-100	50-100	50-100	1X or 2X	☺
Hin6I (HinPI)	G↓CGC	Tango	50-100	50-100	20-50	20-50	100	50-100	1X or 2X	☺
HincII (HindII)	GTY↓RAC	Tango	50-100	50-100	20-50	50-100	100	50-100	1X or 2X	☺
HindIII	A↓GAGCTT	R	0-20	20-50	0-20	100	50-100	50-100	1X or 2X	☺
HinfI	G↓AATC	R	0-20	20-50	50-100	100	50-100	50-100	1X or 2X	☺
HpaI	C↓GGG	Tango	50-100	50-100	0-20	20-50	100	20-50	1X or 2X	☺
HphI	GGTGA(8/7)↓	B	100	0-20	0-20	0-20	20-50	0-20	1X	☺
Hpy8I (MjaI)	GTN↓NAC	Tango	50-100	50-100	0-20	20-50	100	50-100	1X or 2X	☺
HpyF3I (DdeI)	C↓TNAG	Tango	20-50	20-50	20-50	20-50	100	50-100	1X or 2X	☺
HpyF10VI (MwoI)	GCNNNN↓NNGC	Tango	0-20	0-20	0-20	0-20	100	50-100	1X or 2X	☺
KpnI	GTAC↓C	Unique	20-50	0-20	0-20	0-20	20-50	0-20	1X	☺
Kpn2I (BspEI) (55°C)	T↓CCGA	Tango	50-100	50-100	0-20	20-50	100	50-100	1X or 2X	☺
KspAI (HpaI)	GT↓AAC	B	100	50-100*	20-50	20-50	100*	50-100	1X* or 2X*	☺
LglI (SapI)	GCCTC(11/4)↓	Tango	20-50	50-100	20-50	20-50	100	20-50	1X or 2X	☺
Lsp1109I (BbvI)	GCAGC(8/12)↓	Unique	0-20	20-50*	50-100*	100*	20-50*	20-50*	1X* or 2X*	☺
LweI (SnaNI)	GCATC(5/9)↓	Tango	0-20	0-20	0-20	20-50	100	20-50	1X or 2X	☺
MauBI	CG↓CGGCG	Tango	0-20	0-20	0-20	0-20	100	0-20	1X	☺
MbiI (BsrBI)	CGGTC(3/3)↓									