

**Thermo Scientific
Molecular BioProducts Pipette Tips**



find the proper tip
for your pipette

find the proper tip

for your pipette

Offering a broad selection of pipette tips to fit virtually any pipette in your lab, Thermo Scientific MBP pipette tips provide precise sample delivery. From raw materials to the finished product, all components of each manufacturing lot are fully traceable and produced in a controlled environment. Every tip delivers **proven quality** and is certified free of RNase, DNase, DNA, ATP and endotoxin contamination.

These tips are specifically designed and quality tested to meet manufacturing specifications for optimum results when used with Gilson, Eppendorf, Sartorius Biohit and Mettler Toledo Rainin as well as other pipette brands.

MBP® pipette tips deliver:

Environmentally Responsible Packaging

Offering a wide variety of packaging options ranging from racked, individually wrapped, bulk and reload.

Our reload system uses recycled components and 30% less plastic!



Low Retention Technology

The hydrophobic inner surface significantly reduces sample loss due to adhesion and maximizes sample recovery.

Pipette Tip without Low Retention Technology




MBP Pipette Tip with Low Retention Technology



Ensured Protection

ART® self-sealing barrier tips prevent cross contamination













NEW!

Download the **Pipette Tip Selector Apps for iPhone®/iPad®, Android® Today!**

For immediate access to the mobile app of your choice, visit


www.thermoscientific.com/tipselector





Item Number	Description	Packaging	Length (mm)
10 uL			
			
● 3502	10, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 10 trays/pack	31.3
3501	10, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 10 trays/pack	31.3
3504	10, Non-Filtered Tip	Stack Rack, Non-sterile, 96 tips/tray, 2 stacks of 5 trays/pack	31.3
3500	10, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	31.3
● 2139-HR	10, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	31.3
2139-RI	10, ART Barrier Tip	Reload, Filtered, Sterile, 96 tips/tray, 10 trays/pack	31.3
2139B	10, ART Barrier Tip	Bulk, Filtered, Non-sterile, 1000 tips/bag	31.3
20 uL			
			
● 3522	20, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 10 trays/pack	45.7
3521	20, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 10 trays/pack	45.7
3520	20, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	45.7
● 2149P-HR	20, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	50.7
2149P-RI	20, ART Barrier Tip	Reload, Filtered, Sterile, 96 tips/tray, 10 trays/pack	50.7
2149PB	20, ART Barrier Tip	Bulk, Filtered, Non-sterile, 1000 tips/bag	50.7
200uL			
			
● 3552	200, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 10 trays/pack	50.7
3551	200, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 10 trays/pack	50.7
3554	200, Non-Filtered Tip	Stack Rack, Non-sterile, 96 tips/tray, 2 stacks of 5 trays/pack	50.7
3550	200, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	50.7
● 2069-HR	200, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	58.8
2069-RI	200, ART Barrier Tip	Reload, Filtered, Sterile, 96 tips/tray, 10 trays/pack	58.8
2069B	200, ART Barrier Tip	Bulk, Filtered, Non-sterile, 1000 tips/bag	58.8
300uL			
			
● 3772	300, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 10 trays/pack	52.0
3771	300, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 10 trays/pack	52.0
● 2070	300, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	58.4
2070B	300, ART Barrier Tip	Bulk, Non-sterile, 1000 tips/bag	58.4
1000uL			
			
● 3982	1000, Non-Filtered Tip	Racked, Non-sterile, 100 tips/tray, 8 trays/pack	69.9
3981	1000, Non-Filtered Tip	Racked, Sterile, 100 tips/tray, 8 trays/pack	69.9
3984	1000, Non-Filtered Tip	Stack Rack, Non-sterile, 100 tips/tray, 2 stacks of 5 trays/pack	69.9
3980	1000, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	69.9
● 2179-HR	1000, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	88.9
2179-RI	1000, ART Barrier Tip	Reload, Filtered, Sterile, 96 tips/tray, 10 trays/pack	88.9

● denotes tip being shown


For use with Gilson, Eppendorf, Sartorius Biohit and Mettler Toledo Rainin as well as other pipettes brands.


Item Number	Description	Packaging	Length (mm)
Extended Length Tips			
			
● 3512	10 REACH, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 10 trays/pack	43.4
3511	10 REACH, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 10 trays/pack	43.4
3510	10 REACH, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	43.4
● 2140-HR	10 REACH, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	43.3
2140-RI	10 REACH, ART Barrier Tip	Reload, Filtered, Sterile, 96 tips/tray, 10 trays/pack	43.3
2140B	10 REACH, ART Barrier Tip	Bulk, Filtered, Non-sterile, 1000 tips/bag	43.3


			
● 3542	200 REACH, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 8 trays/pack	90.9
3541	200 REACH, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 8 trays/pack	90.9
3540	200 REACH, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	90.9
● 2160P	200 REACH, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 8 trays/pack	90.9


			
● 3792	1000 REACH, Non-Filtered Tip	Racked, Non-sterile, 100 tips/tray, 8 trays/pack	118.5
3791	1000 REACH, Non-Filtered Tip	Racked, Sterile, 100 tips/tray, 8 trays/pack	118.5
● 2079	1000 REACH, ART Barrier Tip	Racked, Filtered, Sterile, 100 tips/tray, 8 trays/pack	118.5

Specialty Tips

			
● 3651	10 GEL, Non-Filtered Tip	Racked, Non-sterile, 200 tips/tray, 1 trays/pack	76.7
● 2155P	20 GEL, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	68.6

			
● 3621	200 GEL, Non-Filtered Tip	Racked, Non-sterile, 200 tips/tray, 1 trays/pack	82.6
● 2155	100 GEL, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 10 trays/pack	68.6

			
● 3702	200, Wide Bore, Non-Filtered Tip	Racked, Non-sterile, 96 tips/tray, 8 trays/pack	91.0
3701	200, Wide Bore, Non-Filtered Tip	Racked, Sterile, 96 tips/tray, 8 trays/pack	91.0
3700	200, Wide Bore, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	91.0
● 2160G	200, Wide Bore, ART Barrier Tip	Racked, Filtered, Sterile, 96 tips/tray, 8 trays/pack	91.0

			
● 3592	1000, Wide Bore, Non-Filtered Tip	Racked, Non-sterile, 100 tips/tray, 8 trays/pack	85.4
3591	1000, Wide Bore, Non-Filtered Tip	Racked, Sterile, 100 tips/tray, 8 trays/pack	85.4
3590	1000, Wide Bore, Non-Filtered Tip	Bulk, Non-sterile, 1000 tips/bag	85.4
● 2079G	1000, Wide Bore, ART Barrier Tip	Racked, Filtered, Sterile, 100 tips/tray, 8 trays/pack	85.4

● denotes tip being shown



Download the **NEW Pipette Tip Selector Apps** for iPhone®/iPad®, Android® Today!

For immediate access to the mobile app of your choice visit www.thermoscientific.com/tipselector

© 2012 Thermo Fisher Scientific Inc. All rights reserved. iPad and iPhone are trademarks of Apple Inc. Android is a trademark of Google Inc. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details

North America: USA/Canada 800 345 0206 or 800 522 7763 3
 Europe: Austria +43 1 801 40 0, Belgium +32 53 73 42 41, France +33 2 2803 2180, Germany national toll free 0800 1 536 376, Germany international +49 6184 90 6000, Italy +39 02 02 95059 552, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (495) 739 76 41, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203
 Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,
 Other Asian countries +852 2885 4613, countries not listed: +49 6184 90 6940 or +33 2 2803 2180

Thermo
 SCIENTIFIC
 Part of Thermo Fisher Scientific

BRHPTMBP0912